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EN English

Keep out of reach of children. Read instructions before use.

Intended Use

The device is intended for cosmetic skin care. Microneedling supports skin renewal and enhances the absorption of cosmetic products. For cosmetic use only.

Instructions for Use

Thoroughly cleanse and disinfect the skin. Fill the reservoir (10 ml) with your chosen serum. Set the desired needle depth using the adjustable dial (0.25 / 0.50 / 1.0 / 1.5 mm). Gently press the stamp perpendicular to the skin without excessive pressure. Move 4–5 times in each direction: horizontally, vertically, and diagonally. Maintain a uniform speed and pressure. Apply a soothing product after treatment.

Recommended Depth

0.25 mm – sensitive areas
 0.50 mm – forehead, maximum depth for home use
 1.0 mm – under professional supervision only
 1.5 mm – under professional supervision only
 Depths exceeding 0.5 mm are not intended for home use.

Safety Warnings

Consult a doctor or dermatologist before first use. Do not use: on injured or irritated skin in case of skin infection or inflammation on active acne if you have blood clotting disorders

during pregnancy Do not share the device with another person.

Possible Reactions

Redness or sensitivity may occur. These effects usually subside within 48 hours. If a persistent or severe reaction occurs, consult a doctor.

Hygiene

The device is intended for personal single use. For hygienic reasons, we recommend disposing of the stamp after use. Do not share the device with another person.

Storage

Store in a dry place and keep out of reach of children.

Disposal

Dispose of the device as sharps waste in accordance with local regulations.

Notice

This product is not intended to treat medical conditions.

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